() R1 () R2 () R3 (X) PG0 () PG1 () Estagiário () Tecnólogo () PIBIC Last Name - Fasolo First Name - Leonardo Middle - Reichmann

Service (sector) Glaucoma Nº CEP

## Sensitivity and specificity of CSLO, SLP, FDT and SWAP in primary open-angle glaucoma: univariate and multivariate analyses.

L.R. Fasolo1, I.M. Tavares1,2, R.B. Fernandes<sup>1</sup>, L.A. S. Melo Jr1, A. Barbosa1, A.Paranhos Jr1, J.A. Prata Jr1,3.

1Vision Institute, Federal University of Sao Paulo, Sao Paulo, Brazil; 2Federal University of Sergipe, Aracaju, Brazil; 3Faculdade de Medicina do Triângulo Mineiro, Uberaba, Brazil.

Purpose: To analyze the diagnostic value of confocal scanning laser ophthalmoscopy (CSLO / Heidelberg Retina Tomography – HRT II), frequency doubling technology perimetry (FDT), scanning laser polarimetry (SLP - GDx) and short-wavelength automated perimetry (SWAP) in normal and glaucomatous patients. Methods: Sixty glaucomatous patients and 60 normal controls were studied. Sensitivity and specificity were analyzed by binary logistic regression. Results: Sensitivity and specificity were, respectively, 55.0% and 98.3% for HRT, 81.6% and 90.0% for FDT, 71.7% and 96.6% for SWAP and 71.7% and 68.3% for GDx. Multiple logistic regression showed statistical significance for HRT (p = 0.009), FDT (p =(0.002) and SWAP (p = 0.003). GDx was not statistically significant (p = 0.098). The probability of positive glaucoma diagnosis was 96.99% when HRT and FDT were abnormal, 99.0% when HRT and SWAP were abnormal, and 99.0% when FDT and SWAP were abnormal. The probability of having glaucoma was 99.9% when all exams were outside normal limits, 76.6%, 44.1% and 70.4% when HRT, FDT and SWAP were outside normal limits, respectively. Conclusion: Diagnostic indices varied among tests. FDT showed higher sensitivity, and HRT and SWAP higher specificity. Better indices were obtained when two or more tests were considered together. Probability of positive glaucoma diagnosis was almost 100% were HRT, FDT and SWAP, all together, outside normal limits.